

Instructions for Operation

The BCX Ultra Deluxe may be used as a Plasma Carrier Modulator alone, as a TENS Stimulator alone, controlling a number of other Ultra Accessories, or a combination of Accessories.

Differential Plasma Carrier Modulation requires BCX Ultra Deluxe Ray tubes

TENS Stimulation requires Ultra Electrodes, such as Hand Cylinders, Foot Plates, or Sticky Patches

External Accessories plug into the BCX Ultra

Running a standard Program:

1. Plug the power supply for the BCX Ultra Deluxe into the rear of the BCX Ultra Deluxe, and the other end into an AC wall outlet.
2. Press On/Off momentarily to Power On
3. Press Change+ or Change- and OK for Prestored or Custom Program
4. Press Change+ or Change- and OK to select Program number
5. Press OK to use Defaults
6. Press Change+ or Change- and OK for Electrodes, Ray tubes, any Accessory, or a Combination
7. Press Change+ or Change- and OK to Select Run Time
8. Plug in Electrodes and Accessory, if used, and press OK
9. Press Start
10. Program starts, displaying Sequences, Frequencies, and remaining time
11. Press Pause to Pause or Continue
12. Auto Shutoff at the end of the program, if selected
13. Press On/Off momentarily to manually shutoff